

S-BIO PrimeSurface® dimensions to support instrument setup

Plate Type	Item number	Well Volume (μL)	Well Depth (mm)	Well Diameter top/bottom (mm)	Plate Length (mm)	Plate Width (mm)	Plate Height (mm)	A1 Row offset (mm)	A1 Column Offset (mm)	Well Center to Well Center spacing (mm)	Flange or Skirt Height (mm)	Stack Height (mm)	Well Bottom Elevation (mm)	Well Bottom Thickness (mm)	Distance to Bottom of Plate (mm)
96 well	MS-9096UZ	300	10	7.0/6.1	127.5	85.9	14	11.45	14.25	9	3.5	12.9	4	1.1	2.9
96 well	MS-9096MZ	201	10	6.7/3.87	127.5	85.5	14	11.25	14.25	9	3.5	13	4	1	3
96 well	MS-9096VZ	298	11.3	6.5/6.17	127.5	85.2	14	11.1	14.25	9	3.5	13	2.7	1	1.7
96 well(Slit)	MS-9096SZ	100	7	6.34/6.17	127.76	85.48	20	11.24	14.38	9	3.5	19	2.7	1	1.7
384 well	MS-9384UZ	100	12	3.3/2.92	127.76	85.48	14.4	8.99	12.13	4.5	2.5	13	2.4	0.8	1.6